


ASTRONOMY 101

Spectra of Gases and Solids

September 26th, 2024



Housekeeping

Lab Rubric, GIMP installation and usage instructions, Spectra images and Submission Portal are on Brightspace

Labs are due October 3rd at 11:59 PM

Housekeeping Continued

Attendance at Labs is Mandatory!

Lab Recordings will be posted on Brightspace for your reference but are not meant as a substitute for class attendance

Late Labs are generally not accepted, please email me if you need an extension for a valid unforeseen reason such as illness or a family emergency

Please Review the Lab Report Expectations in the Lab Manual

Lab Report Overview

Objective/Introduction: One or two sentences for the objective stating the goal(s) of the lab.

The Introduction should be several paragraphs (at least half a page) discussing the scientific concepts of the lab, including any relevant historical context may be included but should not comprise most of the introduction

Procedure: Short paragraph walking through the procedure or steps of the lab.

Observations/Tables/Graphs: All tables, graphs and sketches will go in this section.

Calculations: Sample calculations for each equation used will go in this section

Answers: Answer the the questions in the lab manual in this section.

Lab Report Overview Continued


Discussion:

Arguably the most important section of the lab report, so please do not skip this section!
The discussion should be several paragraphs long (unless otherwise stated) discussing your results and summarizing the following

Did you meet the objective?

- **Why or why not?**
- **Were any assumptions made that may affect your results?**
- **What are the sources of error if any?**
- **What is the broader context?**

Conclusions/References: A short paragraph restating your results where you cite any references used in your lab report (You're expected to use outside sources for each lab). Please also include in text citations as well.



Spectra Lab

Learning Objectives

- Nature of light → electromagnetic radiation
- Generation of a spectrum
- Types of spectra (continuum, absorption and emission)
- Measurement of wavelengths of atomic spectral lines
- Identification of an unknown element by its spectrum




Spectroscopy - the study of the interaction between matter and electromagnetic radiation

From spectral lines astronomers can determine:


- the presence and density of an element
- temperature from the wavelength distribution of the continuum (Wein's Law)
 - add luminosity to get the mass and size
- internal temperature from spectral line width
- the magnetic field from line splitting and polarization
- stellar winds from x-ray emissions
- orbital period and velocity from line movement over time (Doppler Shift)
- the physical changes in the star from line intensity changes over time
- the material around stars - the interstellar medium (ISM).

Spectroscopy is a fundamental tool used to study the universe.


Hydrogen Emission Spectrum



Light Refraction Through a Prism




Prisms



© 2011 Encyclopædia Britannica, Inc.


In 1666, Sir Isaac Newton passed a beam of sunlight through a glass prism producing the visible colour spectrum on a wall.



He then projected light through a second prism to create white light.

Conclusion: the colour spectrum projected is the result of the prism interacting with already-coloured light rather than prism itself generating the colour.

Electromagnetic Waves




A Few Light Facts

Light waves:


- are transverse electromagnetic waves
- carry energy
- are emitted and absorbed by matter
- travel at 300 000 km/s in a vacuum and slower in a more dense medium
- can be reflected and refracted
- travel through a vacuum and some matter

The Electromagnetic Spectrum




velocity = frequency x wavelength = 299 792.458 km/s (speed of light in a vacuum)
~ 300 000 km/s

The Visible Spectrum




The **Ångström** (Å) is a unit of length equal to 10^{-10} m (one ten-billionth of a metre) named after the Swedish physicist Anders Jonas Ångström (1814–1874).

Emission and Absorption Spectra




Continuous spectrum is produced by hot, dense objects

- heated metal, a star, humans, ice cubes and all matter $> 0\text{ K}$
- caused by internal motion and collisions of the particles (*temperature*)




← Increasing Energy








Incandescent Spectrum

Cool tungsten filament
 $T \approx 1000 \text{ K}$

Hot tungsten filament
 $T \approx 2800 \text{ K}$ (melts at 3695 K)




Star Temperature and Colour




Colors are exaggerated


Emission and Absorption Spectra




Emission and Absorption Spectra






Emission and Absorption Spectra



Emission and Absorption Spectra





Gratings and Prisms



PART I: VISUAL OBSERVATION OF DIFFERENT TYPES OF SPECTRA.

- At the front of the lab, there are several different light sources. You will need to sketch & describe the spectra of each one.
- Focus on the prominent lines and their approximate spacing/brightness
- We will allow you to take photos. HOWEVER, there are some things you need to be careful of:
 - You need to make sure you actually photographed the spectrum you're interested in and not one from a nearby source
 - It's quite easy for the sensor on your camera to be saturated by bright sources. In order to get the necessary detail it may be best to sketch
 - The camera may have a hard time picking up any lines present in the solar spectrum

Source 1: Incandescent Bulb Spectra



The horizontal lines in the spectrum are from the bulb supports, they aren't absorption lines.

Take a picture of the spectrum at **low** and **high** power.

Where is it **brightest** (most intense) in each case?

Describe the spectra. What kind of spectra are they?

Source 2: Fluorescent Tube Spectrum




Fluorescent Lights are phosphor-coated tubes filled by low-density Mercury vapour


Mercury vapour releases UV photons → Hit phosphor → Phosphor glows in Visible light

Sketch & Describe the spectrum. What kind of spectrum is it?

Source 3: Sunlight




The Solar Spectrum



Designation	Element	Wavelength (nm)
y	O ₂	898.765
Z	O ₂	822.696
A	O ₂	759.370
B	O ₂	686.719
C	H _α	656.281
a	O ₂	627.661
D1	Na	589.592
D2	Na	588.995
D3 or d	He	587.5618
e	Hg	546.073
E	Fe	527.039
b1	Mg	518.362
b2	Mg	517.270
b3	Fe	516.891
b4	Mg	516.733

Designation	Element	Wavelength (nm)
c	Fe	495.761
F	H _β	486.134
d	Fe	466.814
e	Fe	438.355
G'	Hγ	434.047
G	Fe	430.790
G	Ca	430.774
h	Hδ	410.175
H	Ca ⁺	396.847
K	Ca ⁺	393.366
L	Fe	382.044
N	Fe	358.121
P	Ti ⁺	336.112
T	Fe	302.108
t	Ni	299.444

Sources 4a)-4i): Discharge Tubes



At the front are several discharge tube “carousels.” 3 of them have tubes of:


- a) Air
- b) Argon (Ar)
- c) Carbon Dioxide (CO₂)
- d) Helium (He)
- e) Hydrogen (H)
- f) Neon (Ne)
- g) Nitrogen (N)
- h) Water Vapour (H₂O)

The fourth carousel has an “unknown” gas
Sketch & Describe each source. What kind
of spectrum is it?

Identify the Mystery Element!

PART II: SPECTROGRAPH CALIBRATION USING HELIUM TO ID HYDROGEN


- In this part of the lab, you will be using the known wavelengths of Helium to determine the wavelengths of the first three Hydrogen Balmer lines
- This will be done using **calibrated images** of the two elements
 - The position of a line in a diffraction grating spectrum is dependent on **1) The wavelength** of the light, **2) The properties of the grating**, and **3) the distance to the source**
 - **As long as nothing is changed between the two images** besides the element, we can use the helium image to convert between **line displacement in pixels** and **wavelength in angstroms**



PART II: SPECTROGRAPH CALIBRATION

HOW ASTRONOMERS ID WAVELENGTHS

- Images of spectra are recorded in a digital image → grid of “pixels” of set size
- Measurable quantity: separation
- Line separation from 0th order is proportional to wavelength of light
- Same equipment = same proportional separation for all spectra
- If we take a picture of a known spectrum, eg. Discharge tube, we can use it to map pixels → wavelength
- The lines won't be the same, but the *relationship* between separation and wavelength will → *Calibration*






We will be using the **Measurement tool** in GIMP to measure the displacement of lines in the image



It looks like a compass

If you can't find it, hold
Shift and type "m"
Or
Tools > Measure

Calibration Graph Using Helium Spectral Lines



Calibration Graph Using Helium Spectral Lines

